



MEDI WAVES

Oxygen Concentrator

Model No. MEDAIR
05Ltr-10Ltr



- 
Imported Molecular Sieve
- 
Oil free Compressor
- 
24 hours continuous working
- 
Wise alarm function
- 
High purity oxygen ≥90%
- 
Three-channel oxygen filter
- 
Superior Oxygen atomizing technology
- 
Optional Oxygen purity alarm

- ♦ America PSA technology offer nature oxygen
- ♦ France Imported molecular sieve bed
- ♦ Reliable and durable oil free compressor
- ♦ Available 24 hours continuous working
- ♦ Self diagnostic system with error code indication
- ♦ Optional dual flow and flow splitter for multi people use

Functions

Power off alarm, Overload protection, High/Low pressure alarm, Temperature alarm, Error code indication, Nebulizer, Oxygen purity alarm



SPA Oxygen Concentrator

Oxygen concentrator can not only be for medical use, but also for healthcare and beauty SPA. Konsung SPA oxygen concentrator is specially designed with big oxygen flow and higher pressure. By the spray, it can make your skin breath the fresh high purity oxygen, refresh the cell viability, bring the second wind for your skin.

Specification

Oxygen Purity	Flow Range	Max Noise	Input Voltage	Output Pressure	Power(a)	Weight	Size(mm)	Dimension for outer packing(mm)
93±3%	0-10L	≤52dB	220V/110V	30-70kPa	≤360VA	18kgs	380X320X590	465X375X655



MEDI WAVES INC.

H.O. : B-68, G.T. Karnal Road, Industrial Area, Delhi-110033 (India)

Works : Plot No. 92, Sector-57, Phase IV, HSIIDC Kundli, Sonapat, Haryana-131001 (India) Ph. : 0130-6542024

Website : www.medi-waves.com