


Elite Series


 Easy to maintain with modular design

 Easy replacement of filters

 Neat and innovative design

 Oxygen percentage display

 Timer and time accumulation

 Audible and visual alarm:

- Low oxygen
- Power failure
- Compressor failure
- Low & high flow



Specifications

Model	DO2-5E Series
Oxygen Flow	1-5 LPM
Power Consumption	300W
Output Pressure	42 dB(A)
Net Weight	15 Kgs
Operating Temp.	5°C- 40°C
Dimension/ Packing Size	38*26*58 cm/ 43*32*63 cm
Equip. Classification	Class II Type B
Power Supply	230V, 50Hz or 115V, 60Hz
Optional Function	Oxygen percentage display & low oxygen alarm Low & high flow alarm Nebulizer