

MEDICAL OXYGEN CONCENTRATOR

Model: DO2-5A, DO2-8A, DO2-10A

>> Brief Introduction:

They require very few maintenance among their lifetime. All inner components are modularization designed, very easy to be maintained by general persons. The work performance is very steadily.



Model No	DO2 - 5A	DO2 - 8A	DO2 - 10A
Oxygen flow:	1--5 LPM	1--8 LPM	1--10LPM
Oxygen purity:	93%±3%		
Output pressure	0.05 - 0.07Mpa		
Power Supply:	AC 220V±22V / 50Hz or 110V ± 10 V /60Hz		
Operation noise(dB):	< 50dB (A)		
Power consumption:	320W	540W	720W
Net weight:	15.9 kg	19.2 kg	19.8 kg
LCD display	Normal, total running time, high/low pressure, low oxygen (optional), service required		
Alarm	Power failure, high/low pressure		
Operating temp:	5°C--40°C		
Dimensions (L × W × H):	43 X 32 X 62 cm		
Equip. Classification:	Class II Type B		
Optional	1- Nebulizer 2- Splitter (2-way, 3-way, 4-way, 5-way) 3- Low oxygen alarm and display 4- Double flow outlets		



Single Flow



Double Flow

>> Highlights

- ① Lithium molecular sieve.
- ② GVS imported filters
- ③ USA imported master board.
- ④ Modular design
- ⑤ Packing box with 10mm thickness