

## SPECIFICATIONS

1. Power supply: AC230V, 50Hz/ AC220V, 60Hz/AC120V, 60Hz; Power: 390VA
2. Sound level:  $\leq 46\text{dB (A)}$
3. Maximum recommended flow: 5L/min
4. Flow Range at Outlet Pressure: 0.5~5L/min
5. Oxygen Concentration: 0.5~5L/min ,93% $\pm$ 3%
6. Output Pressure : 38kPa $\pm$ 5kPa;
7. Release Pressure by machine operation :250kPa $\pm$ 50kPa;
8. Net Weight:16.0kg, Gross weight: 18.5kg;
9. Product Dimension: 33 $\times$ 26 $\times$ 54 cm, Package Dimension: 37 $\times$ 43 $\times$ 70 cm;
10. Warranty:  
Long-life compressor designed for 10000 hours of maintenance-free use

## SAFETY SYSTEM:

1. Current overload or line surge shutdown.
2. High temperature compressor shutdown.
3. High pressure alarm shutdown.
4. Low pressure alarm shutdown.
5. Low Oxygen Concentration alarm.

Lead Time  
**0-10**  
DAYS



Model# JM-5A Ni

