

CARER+

New upgrade Intelligent oxygen

P Series Medical Oxygen Concentrator



P5W
New
upgrade



Medical Oxygen Concentrator



CR-P5W

- **Super HD large screen**
Intelligent adjustment of brightness
- **Low noise operation**
Have a quiet time inhaling oxygen
- **Filter replacement**
Intelligent reminder
- **One camera dual use**
Oxygen/Atomization
- **Voice operation prompt**

Power supply AC220V±22V 50Hz±1Hz

Power 450VA

Noise 45dB(A)*

Net weight 16.6kg

Oxygen flow rate 0.5-5L/min

Oxygen concentration 93%±3%

Dimensions 33x34x53cm

Package size 42x40x61cm

*CARER's lab test data may vary from product to product.