

DIGITAL BABY WEIGHING SCALE



Specification

MODEL NO.	P-016
MEASURING CAPACITY	20 Kgs.
GRADUATION	05 gm
CALIBRATION	AUTO
ZERO TRACKING	AUTO
HOLD FACILITY	AVAILABLE
MEASURING PAN	Large Weighing Pan of ABS/SS of size 55x25cms (aprox.)
POWER SUPPLY	Works on Main with 220V (AC)
BACKUP	In-Build Battery backup of 4 Hrs.



MEDI WAVES INC.

ISO 9001:2015 CERTIFIED

WHO-GMP CERTIFIED | ISO 14001:2004 CERTIFIED | OHSAS 18001:2007

ISO 13485 : 2016 CERTIFIED | CE CERTIFIED

H.O. : B-68, G.T. Karnal Road, Industrial Area, Delhi-110033 (India)

Works : Plot No. 92, Sector-57, Phase IV, HSIIDC Kundli,
Sonepat, Haryana-131001 (India) Ph. : 0130-6542024

Helpline No. : 011-42384365

E-mail : sales@medi-waves.com | anandmedical@gmail.com

Website : www.medi-waves.com