



MEDI-WAVES INC.

An ISO 9001 : 9002 certified Co.

**B-68, G T KARNAL ROAD INDUSTRIAL AREA, DELHI.110033 INDIA,
Email: Sales@medi-waves.com <http://www.medi-waves.com>
Tel:91-11-42384227,42384365, Fax: 91-11-42384227**

Electric Suction Apparatus



- Shock Proof, Light Weight Fiber Glass Chassis.
- Noise level of suction apparatus is 50 dB +_ 03 dB.
- Ideal for medical & surgical procedures.
- Rotary Vane or Diaphragm Vacuum Pump with -600 + -10mm Hg.
- Free air displacement 18 ~ 20 ltrs/min.
- Electrical Requirement- 220~230VAC, 50 Hz, 1 Phase.
- 0.125 H.P. Electric Motor.
- 50 mm Vacuum Gauge.
- Non collapsible PVC tubing.
- 2 X 1.0 Ltrs. Glass jars with overflow safety.
- Bacterial filter fitted on back.