

MEDI WAVES INC.

for a healthhier tomorrow...

FETAL DOPPLER

FD-200



FEATURES

- LCD / LED DISPLAY (45MM X 22MM OR MORE)
- IT HAS ULTRASONIC EMITTING FREQUENCY: 2 – 2.5 MHZ AND POWER: < 5M WATT /CM2
- OVERALL SENSITIVITY AT THE DISTANCES 120 – 200 MM FROM THE FACE OF THE PROBE: 90DB
- AUTOMATIC FHR TEST RANGE 50 – 210 BPM OR WIDER (BEAT / MINUTE)
- RESOLUTION: 1 BPM
- AUDIO OUTPUT : 0.5 – 1 WATT
- WITH ELECTRIC CHARGER & RECHARGEABLE BATTERY (LI – ION/ALKALINE)
- GOOD QUALITY WATER PROOF CASE FOR CARRYING PURPOSE
- GOOD QUALITY WATER PROOF PROBE & GEL



MEDI WAVES INC.



ISO 900 : 2015 CERTIFIED
WHO-GMP CERTIFIED | ISO 14001 : 2004 CERTIFIED | OHSAS 18001 : 2007
ISO 13458 : 2016 CERTIFIED | CE CERTIFIED

H.O. : B-68, G.T. KARNAL ROAD, INDUSTRIAL AREA, DELHI - 110033 (INDIA)

Works : Plot No. 92, Sector-57, Phase IV, HSIIDC Kundli,
Sonapat, Haryana-131001 (India) Ph. : 0130-6542024

Helpline No. : 011-42384365

E-mail : sales@medi-waves.com | anandmedical@gmail.com

Website : www.medi-waves.com