



FEATURES

- Durable length board, accurate and safe for infant or child being measured.
- Infant length board Specifications (Infantometer).
- A firm, flat horizontal surface.
- Measuring tape will have 0.1 cm (1mm) increments.
- Fixed tape is stable and easily read.
- An immovable headpiece at right angle to the tape, A movable foot piece, perpendicular to tape, length board have a measuring range to at least 40 to 99 cm.

Specification

MODEL NAME	P-017(A)
UNIT OF MEASURING	cm
MEASURING RANGE	40 -99 cm
GRADUATION	1mm



MEDI WAVES INC.

(An ISO 13485 : 2003 Certified Co.)

H.O. : B-68, G.T. Karnal Road, Industrial Area, Delhi-110033 (India)

Works : Plot No. 92, Sector-57, Phase IV, HSIIDC Kundli,
Sonapat, Haryana-131001 (India)

Helpline No. : 011-42384365

Website : www.medi-waves.com

