

Portable Oxygen Concentrator

- Suitable for ambulance
- Oxygen Concentration: 90% - 98%
- Monitor oxygen purity and pressure sensor which capable of real-time monitoring of user activity and sending instant alerts
- Durable and light weight: 2.5 kgs
- Oxygen output: up to a minimum of 2000mlr altitude
- Breath detection sensor and Pulse dose concentrators adjust automatically breath rate
- Outlet Pressure: 45-80KPa
- Pulse Flow rate: Setting:1-5 LPM
- Altering the user of fault conditions: high/lowpressure, low/high flow oxygen, high temperature, power failure, low battery, and no flow
- Sound Level: ≤ 50 dB
- Operated on mains as well as on battery
- Operating Temperature: 0°C to $+ 40^{\circ}\text{C}$
- Relative humidity: $< 85\%$
- Dual battery & Operate for min 10 hrs

MEDIWAVES INC.

Regd. Office :

B-68, G.T. Karnal Road Industrial Area,
Delhi-110033 Ph.: 011-42984227
Web: www.medi-waves.com

MEDAIR - 5 LPM



2.5KG

