



Slow Suction Machine - P-023

Portable suction units is a newly designed oil-free lubrication pump based on extensive design experience and development trend of domestic and overseas products. It is specially designed for absorbing thick liquid such as blood, phlegm etc. It is not only used as aspirator in operation room. But also at home as healthcare product.

Specifications

Features

High negative pressure, large flux, low noise.

It adopts advanced overflow-protection design, which can prevent pus and mucus from leaking into the machine.

Compact, lightweight and easy to carry especially suitable for emergency and outside use.

Voltage : AC220 / 50Hz

Max negative pressure value : more than 0.08
MPA or 550mmHg

Pumping rate : more that 20L per minutes

Reservoir volume : 1000ml

Noise : Less than 65db(A)

Power : 90VA
